

## DOSY (Diffusion Ordered Spectroscopy)

### Background

In DOSY, the spins are gradient-encoded by their z-position in the NMR tube. After a brief diffusion time ( $\Delta$ ) the spins are decoded. If a molecule remains within the same slice of the NMR tube, it will appear at full intensity, but if it diffuses to another slice, its signal will be attenuated as a function of the distance it has diffused. A simplified pulse sequence is shown below:



The peak attenuation is controlled by three variables. Typically, the gradient power ( $G$ ) is varied over the course of the 2D experiment. The diffusion delay ( $\Delta$ ) and gradient length ( $\delta$ ), however, must be optimized.

### Reducing the Spectral Window

1. Insert a sample, lock, shim, and record a normal  $^1\text{H}$  spectrum.
2. Type IEXPNO to move to a new EXPNO. Expand the spectrum so that all the peaks are on the screen, and click the  $\downarrow$  icon to adjust SW and O1. Reduce TD to 16K and SI to 8K.
3. Acquire another spectrum to check that the parameters are good.

### Optimizing the DOSY Parameters:

4. In the ASED window, set PULPRG = StebPgp1s1D, GPz6 = 2 (the percentage of the full gradient to use for this spectrum), GPz7 = -17.13 (homospoil), D20 = 0.1 s ( $\Delta$ -value), P30 = 1800  $\mu\text{s}$  ( $\delta$ -value). Turn the spinner off ("ro off").
5. Obtain a spectrum (rga / zg / efp / apk).
6. Move to a new expno (IEXPNO), change GPz6 to 95 percent of the full gradient, and acquire another spectrum.
7. Compare the two spectra: MDISP then RE  $n$ , where  $n$  is the expno from step 5 (eg, "re 3"). If the intensity of the second spectrum is roughly 1/50 that of the first, all is good. If not, adjust D20 ( $\Delta$ ) or P30 ( $\delta$ ) until they are.

### Running the DOSY

8. Change the pulse program to StebPgp1s. Click the  $\downarrow$  icon under the Pulse Sequence tab to change the data-set to 2D. Increase NS and DS (typically 32 + 4). Set TD(F2,F1) = (16k,  $n$ ), where  $n$  is the number of F1 slices desired (ie, number of G values to use, typically 16). Set SI(F2,F1) to (8k,  $n$ ), and FnMOD = QF.
9. Type DOSY, which will ask about additional parameters. Set the Gradient Ramp using the minimum and maximum values you used in steps 4 and 6 (2% and 95%). Set the number of points to the same  $n$  value used above (16). A linear ramp is most common. Press OK and step back.

### Processing the DOSY Data

(see below if you want numerical diffusion coefficients, rather than a 2-D spectrum)

10. Set PhMod(F2,F1) to (pk, no). Type SetDiffParam to transfer key parameters, and XF2 to transform the data in F2.
11. If necessary, phase the data (check by viewing with the  $\wedge$  button).

13. Process the frequency *vs* time data using DOSY2D, to obtain a frequency *vs* diffusion-coefficient plot.
14. View slices across F1 through the “Slice” / “interactive” command: click the  icon. The slices may be saved to a procno by right-clicking and choosing “extract”.

#### Obtaining Numerical Diffusion Coefficients

The data may be explicitly regressed and printed as a report, rather than graphically. This follows the same principle used for obtaining T1 and T2 values.

- 10'. Set PhMod(F2,F1) to (pk, no).
- 11'. XF2 to transform the data.
- 12'. From the menus, choose Analysis : T1/T2 Relaxation. Follow the flowchart down the left-column of buttons:
  - a. Extract a slice from your data to phase. Choose “spectrum” and use the first slice (default). Phase the data (-button on top row), then save-and-return ( button within phase window).
  - b. Use the flow-chart’s Peaks/Ranges button to define the ranges of the spectrum to regress. Integrate the spectrum (-button) by dragging across the peaks, then export the data back to the T1/T2 routine (-button; export to Relaxation Module).
  - c. Use the flow-chart’s Relaxation Window button to regress the peak intensity *vs*  $\tau$  data. On Fitting Functions step, ( icon), set the Fitting Function to VARGRAD and the List File name to DIFFLIST. The + and – buttons will allow you to change which peak the plot displays data for. The -button will allow you to delete any bad data-points; if you do, use the  and  buttons (left edge of bar) to recalculate the fit.